





## New DOUBLE ROW Foot Plate

- Quicker and easier frame assembly
- Stronger constructs for open foot supports
- Greater freedom for wire attachment

TrueLok is Simple  
TrueLok is Stable  
TrueLok is Versatile





Foot plate ends allow attachment of full rings, 1/2 rings and standard foot plates

Single plate construction provides more options for wire placement

Double row of holes provide more options for wire attachment

M/L reference marks allow for quick frame alignment

Ergonomic design reduces tripping hazards

Posterior reference mark allows for quick frame alignment

Flat posterior of plate supports neutral limb alignment in a supine position

**Available in Popular Sizes:**

Inner Diameter	Part #
160mm	56-13635
180mm	56-13655
200mm	56-13665

Adult & Pediatric Deformity Correction | Foot & Ankle | Bone Growth Stimulation | Fixation | Biologics